

**Table of contents**

\USER

COCOAN

fMRI MBEPi

MPC

[heat\\_1\\_2.7iso\\_ap\\_64ch\\_mb8](#)

\*

\\USER\COCOAN\fmri MBEP\MPC\heat\_1\_2.7iso\_ap\_64ch\_mb8 \*

TA: 6:20 PM: FIX Voxel size: 2.7×2.7×2.7 mmPAT: Off Rel. SNR: 1.00 : efpid

**Properties**

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	On
Start measurements	Single measurement

**Routine**

Slice group	1
Slices	56
Dist. factor	0 %
Position	R6.2 P4.9 F19.8 mm
Orientation	T > C-19.8 > S-1.4
Phase enc. dir.	A >> P
AutoAlign	---
Phase oversampling	0 %
FoV read	220 mm
FoV phase	100.0 %
Slice thickness	2.70 mm
TR	460 ms
TE	27.20 ms
Multi-band accel. factor	8
Filter	Raw filter, Prescan Normalize
Coil elements	HC1-7;NC1

**Contrast - Common**

TR	460 ms
TE	27.20 ms
MTC	Off
Magn. preparation	None
Flip angle	44 deg
Fat suppr.	Fat sat.

**Contrast - Dynamic**

Averaging mode	Long term
Reconstruction	Magnitude
Measurements	812
Delay in TR	0 ms
Multiple series	Off

**Resolution - Common**

FoV read	220 mm
FoV phase	100.0 %
Slice thickness	2.70 mm
Base resolution	82
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off

**Resolution - iPAT**

PAT mode	None
----------	------

**Resolution - Filter Image**

Distortion Corr.	Off
Prescan Normalize	On

**Resolution - Filter Rawdata**

Raw filter	On
Elliptical filter	Off
Hamming	Off

**Geometry - Common**

Slice group	1
Slices	56
Dist. factor	0 %
Position	R6.2 P4.9 F19.8 mm
Orientation	T > C-19.8 > S-1.4
Phase enc. dir.	A >> P
FoV read	220 mm
FoV phase	100.0 %
Slice thickness	2.70 mm
TR	460 ms
Multi-slice mode	Interleaved
Series	Interleaved
Multi-band accel. factor	8

**Geometry - AutoAlign**

Slice group	1
Position	R6.2 P4.9 F19.8 mm
Orientation	T > C-19.8 > S-1.4
Phase enc. dir.	A >> P
AutoAlign	---
Initial Position	R6.2 P4.9 F19.8
R	6.2 mm
P	4.9 mm
F	19.8 mm
Initial Rotation	0.00 deg
Initial Orientation	T > C
T > C	-19.8
> S	-1.4

**Geometry - Saturation**

Fat suppr.	Fat sat.
Special sat.	None

**System - Miscellaneous**

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	On - AutoCoilSelect

**System - Adjustments**

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off

**System - Adjustments**

Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

**System - Adjust Volume**

Position	R6.2 P4.9 F19.8 mm
Orientation	T > C-19.8 > S-1.4
Rotation	0.00 deg
A >> P	220 mm
R >> L	220 mm
F >> H	152 mm
Reset	Off

**System - pTx Volumes**

B1 Shim mode	TrueForm
Excitation	Standard

**System - Tx/Rx**

Frequency 1H	123.253601 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

**Physio - Signal1**

1st Signal/Mode	None
TR	460 ms
Multi-band accel. factor	8

**BOLD**

GLM Statistics	Off
Dynamic t-maps	Off
Ignore meas. at start	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	Off
Spatial filter	Off
Measurements	812
Delay in TR	0 ms
Multiple series	Off

**Sequence - Part 1**

Introduction	Off
Contrasts	1
Flow comp.	No
Multi-slice mode	Interleaved
Free echo spacing	Off
Echo spacing	0.49 ms
Bandwidth	3048 Hz/Px

**Sequence - Part 2**

EPI factor	82
Gradient mode	Performance
Excitation	Standard
RF spoiling	Off

**Sequence - Special**

Excite pulse duration	5120 us
Single-band images	On
MB LeakBlock kernel	Off
MB dual kernel	Off
MB RF phase scramble	Off
SENSE1 coil combine	Off
Invert RO/PE polarity	Off
Disable freq. update	Off
Force equal slice timing	Off
Online multi-band recon.	Online
FFT scale factor	0.90
Physio recording	Off
Triggering scheme	Standard